Public Health
England

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
Bibliografia

1. Saito T, Iinuma Y, Takakura S, Nagao M, Matsushima A, Shirano M, et al. Delayed insertion of blood culture bottles into automated continuously monitoring blood culture systems increases the time from blood sample collection to the detection of microorganisms in bacteremic patients. *J Infect Chemother* 2009;15:49-53.
2. van der Velden LB, Vos FJ, Mouton JW, Sturm PD. Clinical impact of preincubation of blood cultures at 37 degrees C. *J Clin Microbiol* 2011;49:275-80.
3. Clinical Pathology Accreditation (UK) Ltd. Standards for the Medical Laboratory. 2010.
4. Janapatla RP, Yan JJ, Chien ML, Chen HM, Wu HM, Wu JJ. Effect of overnight storage of blood culture bottles on bacterial detection time in the BACTEC 9240 blood culture system. *J Microbiol Immunol Infect* 2010;43:126-32.
5. Lemming L, Holt HM, Petersen IS, Ostergaard C, Bruun B. Bactec 9240 blood culture system: to preincubate at 35 degrees C or not? *Clin Microbiol Infect* 2004;10:1089-91.
6. Ronnberg C, Mildh M, Ullberg M, Ozenci V. Transport time for blood culture bottles: underlying factors and its consequences. *Diagn Microbiol Infect Dis* 2013;76:286-90.
7. Kerremans JJ, Verboom P, Stijnen T, Hakkaart-van Roijen L, Goessens W, Verbrugh HA, et al. Rapid identification and antimicrobial susceptibility testing reduce antibiotic use and accelerate pathogen-directed antibiotic use. *J Antimicrob Chemother* 2008;61:428-35.
8. Kerremans JJ, Goessens WH, Verbrugh HA, Vos MC. Accuracy of identification and susceptibility results by direct inoculation of Vitek 2 cards from positive BACTEC cultures. *Eur J Clin Microbiol Infect Dis* 2004;23:892-8.
9. Edelmann A, Pietzcker T, Wellinghausen N. Comparison of direct disk diffusion and standard microtitre broth dilution susceptibility testing of blood culture isolates. *J Med Microbiol* 2007;56:202-7.
10. Dark P, Dunn G, Chadwick P, Young D, Bentley A, Carlson G, et al. The clinical diagnostic accuracy of rapid detection of healthcare-associated bloodstream infection in intensive care using multipathogen real-time PCR technology. *BMJ Open* 2011;1:e000181.
11. Paolucci M, Landini MP, Sambri V. Conventional and molecular techniques for the early diagnosis of bacteraemia. *Int J Antimicrob Agents* 2010;36 Suppl 2:S6-16.
12. Reimer LG, Wilson ML, Weinstein MP. Update on detection of bacteremia and fungemia. *Clin Microbiol Rev* 1997;10:444-65.
13. Masterson KC, McGowan JE, Jr. Detection of positive blood cultures by the Bactec NR660. The clinical importance of five versus seven days of testing. *Am J Clin Pathol* 1988;90:91-4.
14. Hardy DJ, Hulbert BB, Migneault PC. Time to detection of positive BacT/Alert blood cultures and lack of need for routine subculture of 5- to 7-day negative cultures. *J Clin Microbiol* 1992;30:2743-5.
15. Kurtoglu MG, Bozkurt H, Tuncer O, Kesli R, Berktaş M. Distribution, optimum detection time and antimicrobial susceptibility rates of the microorganisms isolated from blood cultures over a 4-year time period in a Turkish university hospital and a review of the international literature. *J Int Med Res* 2008;36:1261-72.
16. Ozkurt Z, Erol S, Tasyaran MA, Kaya A. Detection of *Brucella melitensis* by the BacT/Alert automated system and *Brucella* broth culture. *Clin Microbiol Infect* 2002;8:749-52.

17. Sumerkan B, Gokahmetoglu S, Esel D. Brucella detection in blood: comparison of the BacT/Alert standard aerobic bottle, BacT/Alert FAN aerobic bottle and BacT/Alert enhanced FAN aerobic bottle in simulated blood culture. Clin Microbiol Infect 2001;7:369-72.
18. National Institute for Healthcare and Clinical Excellence. Antibiotics for early-onset neonatal infection: antibiotics for the prevention and treatment of early-onset neonatal infection. 2012.
19. Labarca JA, Wagar EA, Grasmick AE, Kokkinos HM, Bruckner DA. Critical evaluation of 4-week incubation for fungal cultures: is the fourth week useful? J Clin Microbiol 1998;36:3683-5.
20. Kumar Y, Qunibi M, Neal TJ, Yoxall CW. Time to positivity of neonatal blood cultures. Arch Dis Child Fetal Neonatal Ed 2001;85:F182-F186.
21. Stefani S. Diagnostic techniques in bloodstream infections: where are we going? Int J Antimicrob Agents 2009;34 Suppl 4:S9-12.
22. Rowther FB, Rodrigues CS, Deshmukh MS, Kapadia FN, Hegde A, Mehta AP, et al. Prospective comparison of eubacterial PCR and measurement of procalcitonin levels with blood culture for diagnosing septicemia in intensive care unit patients. J Clin Microbiol 2009;47:2964-9.
23. Fenollar F, Raoult D. Molecular diagnosis of bloodstream infections caused by non-cultivable bacteria. Int J Antimicrob Agents 2007;30 Suppl 1:S7-15.
24. Hautala T, Syrjala H, Lehtinen V, Kauma H, Kauppila J, Kujala P, et al. Blood culture Gram stain and clinical categorization based empirical antimicrobial therapy of bloodstream infection. Int J Antimicrob Agents 2005;25:329-33.
25. Deresinski S. Principles of Antibiotic Therapy in Severe Infections: Optimizing the Therapeutic Approach by Use of Laboratory and Clinical Data. Clinical Infectious Diseases 2007;45:S177-S183.
26. Berild D, Mohseni A, Diep LM, Jensenius M, Ringertz SH. Adjustment of antibiotic treatment according to the results of blood cultures leads to decreased antibiotic use and costs. J Antimicrob Chemother 2006;57:326-30.
27. Sprung CL. Definitions of sepsis--have we reached a consensus? Crit Care Med 1991;19:849-51.
28. Weinstein MP, Towns ML, Quartey SM, Mirrett S, Reimer LG, Parmigiani G, et al. The clinical significance of positive blood cultures in the 1990s: a prospective comprehensive evaluation of the microbiology, epidemiology, and outcome of bacteremia and fungemia in adults. Clin Infect Dis 1997;24:584-602.
29. Jumaa PA, Chattopadhyay B. Pseudobacteraemia. J Hosp Infect 1994;27:167-77.
30. Noskin GA, Suriano T, Collins S, Sesler S, Peterson LR. Paenibacillus macerans pseudobacteremia resulting from contaminated blood culture bottles in a neonatal intensive care unit. Am J Infect Control 2001;29:126-9.
31. Dellinger RP, Levy MM, Carlet JM, Bion J, Parker MM, Jaeschke R, et al. Surviving Sepsis Campaign: international guidelines for management of severe sepsis and septic shock: 2008. Intensive Care Med 2008;34:17-60.
32. Al MM, Darouiche RO. Sepsis syndrome, bloodstream infections, and device-related infections. Med Clin North Am 2012;96:1203-23.
33. Annane D, Bellissant E, Cavaillon JM. Septic shock. Lancet 2005;365:63-78.
34. Cohen J. Sepsis. Medicine 2009;37:562-5.

35. Kerremans JJ, van der Bij AK, Goessens W, Verbrugh HA, Vos MC. Immediate incubation of blood cultures outside routine laboratory hours of operation accelerates antibiotic switching. *J Clin Microbiol* 2009;47:3520-3.
36. Daniels R. Surviving the first hours in sepsis: getting the basics right (an intensivist's perspective). *Journal of Antimicrobial Chemotherapy* 2011;66:ii11-ii23.
37. Paolucci M, Landini MP, Sambri V. How can the microbiologist help in diagnosing neonatal sepsis? *Int J Pediatr* 2012;120:39.
38. Edmond K, Zaidi A. New approaches to preventing, diagnosing, and treating neonatal sepsis. *PLoS Med* 2010;7:e1000213.
39. Brook I. The role of anaerobic bacteria in bacteremia. *Anaerobe* 2010;16:183-9.
40. Finelli L, Livengood JR, Saiman L. Surveillance of pharyngeal colonization: detection and control of serious bacterial illness in low birth weight infants. *Pediatr Infect Dis J* 1994;13:854-9.
41. Lee PY, Holliman RE, Davies EG. Surveillance cultures on neonatal intensive care units. *J Hosp Infect* 1995;29:233-6.
42. Guerti K, Devos H, Ieven MM, Mahieu LM. Time to positivity of neonatal blood cultures: fast and furious? *J Med Microbiol* 2011;60:446-53.
43. Ispahani P, Pearson NJ, Greenwood D. An analysis of community and hospital-acquired bacteraemia in a large teaching hospital in the United Kingdom. *Q J Med* 1987;63:427-40.
44. Plummer M, Schoch PE. *Rothia dentocariosa* bacteremia. *Clinical Microbiology Newsletter* 1995;17:22-4.
45. Auzias A, Bollet C, Ayari R, Drancourt M, Raoult D. *Corynebacterium freneyi* bacteremia. *J Clin Microbiol* 2003;41:2777-8.
46. Roberts FJ. Nontyphoidal, nonparatyphoidal salmonella septicemia in adults. *Eur J Clin Microbiol Infect Dis* 1993;12:205-8.
47. Wilks D, Jacobson SK, Lever AM, Farrington M. Fatal melioidosis in a tourist returning from Thailand. *J Infect* 1994;29:87-90.
48. Collazos J, de Miguel J, Ayarza R. *Moraxella catarrhalis* bacteremic pneumonia in adults: two cases and review of the literature. *Eur J Clin Microbiol Infect Dis* 1992;11:237-40.
49. Wenzel RP. Nosocomial candidemia: risk factors and attributable mortality. *Clin Infect Dis* 1995;20:1531-4.
50. Warnock EW, III, MacMath TL. Primary *Vibrio vulnificus* septicemia. *J Emerg Med* 1993;11:153-6.
51. Khan AM, Albert MJ, Sarker SA, Bhattacharya MK, Azad AK. Septicemia due to *Vibrio cholerae* 0139 Bengal. *Diagn Microbiol Infect Dis* 1995;22:337-8.
52. Bonatti H, Rossboth DW, Nachbaur D, Fille M, Aspöck C, Hend I, et al. A series of infections due to *Capnocytophaga* spp in immunosuppressed and immunocompetent patients. *Clin Microbiol Infect* 2003;9:380-7.
53. Fenner L, Widmer AF, Straub C, Frei R. Is the incidence of anaerobic bacteremia decreasing? Analysis of 114,000 blood cultures over a ten-year period. *J Clin Microbiol* 2008;46:2432-4.

54. Mylotte JM, Tayara A. Blood cultures: clinical aspects and controversies. *Eur J Clin Microbiol Infect Dis* 2000;19:157-63.
55. Riley JA, Heiter BJ, Bourbeau PP. Comparison of recovery of blood culture isolates from two BacT/ALERT FAN aerobic blood culture bottles with recovery from one FAN aerobic bottle and one FAN anaerobic bottle. *J Clin Microbiol* 2003;41:213-7.
56. Paisley JW, Lauer BA. Pediatric blood cultures. *Clin Lab Med* 1994;14:17-30.
57. Jaffe DM. Occult bacteremia in children. [Review] [100 refs]. *Adv Pediatr Infect Dis* 1994;9:237-60.
58. Bouza E, Burillo A, Munoz P. Catheter-related infections: diagnosis and intravascular treatment. *Clin Microbiol Infect* 2002;8:265-74.
59. Hall KK, Lyman JA. Updated review of blood culture contamination. *Clin Microbiol Rev* 2006;19:788-802.
60. Garcia X, Sabatier C, Ferrer R, Fontanals D, Duarte M, Colomina M, et al. Differential time to positivity of blood cultures: a valid method for diagnosing catheter-related bloodstream infections in the intensive care unit. *Med Intensiva* 2012;36:169-76.
61. Mermel LA, Allon M, Bouza E, Craven DE, Flynn P, O'Grady NP, et al. Clinical practice guidelines for the diagnosis and management of intravascular catheter-related infection: 2009 Update by the Infectious Diseases Society of America. *Clin Infect Dis* 2009;49:1-45.
62. O'Grady NP, Alexander M, Burns LA, Dellinger EP, Garland J, Heard SO, et al. Guidelines for the prevention of intravascular catheter-related infections. *Am J Infect Control* 2011;39:S1-34.
63. Kaasch AJ, Rieg S, Hellmich M, Kern WV, Seifert H. Differential time to positivity is not predictive for central line-related *Staphylococcus aureus* bloodstream infection in routine clinical care. *J Infect* 2014;68:58-61.
64. Guembe M, Rodriguez-Creixems M, Sanchez-Carrillo C, Martin-Rabadan P, Bouza E. Differential time to positivity (DTTP) for the diagnosis of catheter-related bloodstream infection: do we need to obtain one or more peripheral vein blood cultures? *Eur J Clin Microbiol Infect Dis* 2012;31:1367-72.
65. Beekmann SE, Diekema DJ, Huskins WC, Herwaldt L, Boyce JM, Sherertz RJ, et al. Diagnosing and reporting of central line-associated bloodstream infections. *Infect Control Hosp Epidemiol* 2012;33:875-82.
66. Janakiraman V. Listeriosis in pregnancy: diagnosis, treatment, and prevention. *Rev Obstet Gynecol* 2008;1:179-85.
67. Mylonakis E, Paliou M, Hohmann EL, Calderwood SB, Wing EJ. Listeriosis during pregnancy: a case series and review of 222 cases. *Medicine (Baltimore)* 2002;81:260-9.
68. Fitzsimmons K, Bamber AI, Smalley HB. Infective endocarditis: changing aetiology of disease. *Br J Biomed Sci* 2010;67:35-41.
69. Moreillon P, Que YA. Infective endocarditis. *Lancet* 2004;363:139-49.
70. Fefer P, Raveh D, Rudensky B, Schlesinger Y, Yinnon AM. Changing epidemiology of infective endocarditis: a retrospective survey of 108 cases, 1990-1999. *Eur J Clin Microbiol Infect Dis* 2002;21:432-7.
71. Brouqui P, Raoult D. Endocarditis due to rare and fastidious bacteria. *Clin Microbiol Rev* 2001;14:177-207.

72. McDonald JR. Acute Infective Endocarditis. *Infect Dis Clin North Am* 2009;23:643-64.
73. Kohler W. The present state of species within the genera *Streptococcus* and *Enterococcus*. *Int J Med Microbiol* 2007;297:133-50.
74. Hogevis H, Alestig K. Fungal endocarditis--a report on seven cases and a brief review. *Infection* 1996;24:17-21.
75. Tacke D, Koehler P, Cornely OA. Fungal endocarditis. *Curr Opin Infect Dis* 2013;26:501-7.
76. Petti CA, Bhally HS, Weinstein MP, Joho K, Wakefield T, Reller LB, et al. Utility of extended blood culture incubation for isolation of *Haemophilus*, *Actinobacillus*, *Cardiobacterium*, *Eikenella*, and *Kingella* organisms: a retrospective multicenter evaluation. *J Clin Microbiol* 2006;44:257-9.
77. Baron EJ, Scott JD, Tompkins LS. Prolonged incubation and extensive subculturing do not increase recovery of clinically significant microorganisms from standard automated blood cultures. *Clin Infect Dis* 2005;41:1677-80.
78. Petti CA, Zaidi AK, Mirrett S, Reller LB. Comparison of Isolator 1.5 and BACTEC NR660 aerobic 6A blood culture systems for detection of fungemia in children. *J Clin Microbiol* 1996;34:1877-9.
79. Holliman R. The immunocompromised patient: infection in cancer and transplantation. *Medicine* 2009;37:522-4.
80. Varani S, Stanzani M, Paolucci M, Melchionda F, Castellani G, Nardi L, et al. Diagnosis of bloodstream infections in immunocompromised patients by real-time PCR. *J Infect* 2009;58:346-51.
81. Algar V, Novelli V. Infections in the immunocompromised host. *Paediatr Child Health* 2007;17:132-6.
82. Beebe JL, Koneman EW. Recovery of uncommon bacteria from blood: association with neoplastic disease. *Clin Microbiol Rev* 1995;8:336-56.
83. Kauffman CA. Diagnosis of histoplasmosis in immunosuppressed patients. *Curr Opin Infect Dis* 2008;21:421-5.
84. Lobmaier IV, Vege A, Gaustad P, Rognum TO. Bacteriological investigation--significance of time lapse after death. *Eur J Clin Microbiol Infect Dis* 2009;28:1191-8.
85. Silver H, Sonnenwirth AC. A practical and efficacious method for obtaining significant postmortem blood cultures. *Am J Clin Pathol* 1969;52:433-7.
86. Pryce JW, Roberts SE, Weber MA, Klein NJ, Ashworth MT, Sebire NJ. Microbiological findings in sudden unexpected death in infancy: comparison of immediate postmortem sampling in casualty departments and at autopsy. *J Clin Pathol* 2011;64:421-5.
87. Wilson SJ, Wilson ML, Reller LB. Diagnostic utility of postmortem blood cultures. *Arch Pathol Lab Med* 1993;117:986-8.
88. Scottish Government. Public Health (Scotland) Act. 2008 (as amended).
89. Scottish Government. Public Health etc. (Scotland) Act 2008. Implementation of Part 2: Notifiable Diseases, Organisms and Health Risk States. 2009.
90. The Welsh Assembly Government. Health Protection Legislation (Wales) Guidance. 2010.
91. Home Office. Public Health Act (Northern Ireland) 1967 Chapter 36. 1967 (as amended).

92. Kanj SS, Kanafani ZA. Current concepts in antimicrobial therapy against resistant gram-negative organisms: extended-spectrum beta-lactamase-producing Enterobacteriaceae, carbapenem-resistant Enterobacteriaceae, and multidrug-resistant *Pseudomonas aeruginosa*. *Mayo Clin Proc* 2011;86:250-9.
93. Ammerlaan HS, Harbarth S, Buiting AG, Crook DW, Fitzpatrick F, Hanberger H, et al. Secular trends in nosocomial bloodstream infections: antibiotic-resistant bacteria increase the total burden of infection. *Clin Infect Dis* 2013;56:798-805.
94. Huttunen R, Syrjanen J, Vuento R, Aittoniemi J. Current concepts in the diagnosis of blood stream infections. Are novel molecular methods useful in clinical practice? *Int J Infect Dis* 2013;17:e934-e938.
95. Rodriguez-Bano J, Navarro MD, Romero L, Muniain MA, de Cueto M, Rios MJ, et al. Bacteremia due to extended-spectrum beta -lactamase-producing *Escherichia coli* in the CTX-M era: a new clinical challenge. *Clin Infect Dis* 2006;43:1407-14.
96. Worth LJ, Slavin MA. Bloodstream infections in haematology: Risks and new challenges for prevention. *Blood Rev* 2009;23:113-22.
97. Vigano EF, Vasconi E, Agrappi C, Clerici P. Use of simulated blood cultures for time to detection comparison between BacT/ALERT and BACTEC 9240 blood culture systems. *Diagn Microbiol Infect Dis* 2002;44:235-40.
98. Rohner P, Pepey B, Auckenthaler R. Comparison of BacT/Alert with Signal blood culture system. *J Clin Microbiol* 1995;33:313-7.
99. Wilson ML, Mirrett S, McDonald LC, Weinstein MP, Fune J, Reller LB. Controlled clinical comparison of bioMerieux VITAL and BACTEC NR-660 blood culture systems for detection of bacteremia and fungemia in adults. *J Clin Microbiol* 1999;37:1709-13.
100. Spanjaard L, Kuijper EJ, Dankert J. Clinical comparison of two commercial blood culture systems. *Eur J Clin Microbiol Infect Dis* 2000;19:881-5.
101. Krisher KK, Gibb P, Corbett S, Church D. Comparison of the BacT/Alert PF pediatric FAN blood culture bottle with the standard pediatric blood culture bottle, the Pedi-BacT. *J Clin Microbiol* 2001;39:2880-3.
102. McDonald LC, Fune J, Gaido LB, Weinstein MP, Reimer LG, Flynn TM, et al. Clinical importance of increased sensitivity of BacT/Alert FAN aerobic and anaerobic blood culture bottles. *J Clin Microbiol* 1996;34:2180-4.
103. Department of Health. Taking blood cultures. 2011.
104. Stohl S, Benenson S, Svirj A, Avidan A, Block C, Sprung CL, et al. Blood cultures at central line insertion in the intensive care unit: comparison with peripheral venipuncture. *J Clin Microbiol* 2011;49:2398-403.
105. Dwivedi S, Bhalla R, Hoover DR, Weinstein MP. Discarding the initial aliquot of blood does not reduce contamination rates in intravenous-catheter-drawn blood cultures. *J Clin Microbiol* 2009;47:2950-1.
106. Byard RW. Arterial blood cultures in disseminated fungal disease. *Pediatr Infect Dis J* 1989;8:728-9.
107. Krumholz HM, Cummings S, York M. Blood culture phlebotomy: switching needles does not prevent contamination. *Ann Intern Med* 1990;113:290-2.
108. Leisure MK, Moore DM, Schwartzman JD, Hayden GF, Donowitz LG. Changing the needle when inoculating blood cultures. A no-benefit and high-risk procedure. *JAMA* 1990;264:2111-2.

109. Spitalnic SJ, Woolard RH, Mermel LA. The significance of changing needles when inoculating blood cultures: a meta-analysis. Clin Infect Dis 1995;21:1103-6.
110. Washington JA, Ilstrup DM. Blood cultures: issues and controversies. Rev Infect Dis 1986;8:792-802.
111. Li J, Plorde JJ, Carlson LG. Effects of volume and periodicity on blood cultures. J Clin Microbiol 1994;32:2829-31.
112. Connell TG, Rele M, Cowley D, Buttery JP, Curtis N. How reliable is a negative blood culture result? Volume of blood submitted for culture in routine practice in a children's hospital. Pediatrics 2007;119:891-6.
113. Mermel LA, Maki DG. Detection of bacteremia in adults: consequences of culturing an inadequate volume of blood. Ann Intern Med 1993;119:270-2.
114. Patel R, Vetter EA, Harmsen WS, Schleck CD, Fadel HJ, Cockerill FR, III. Optimized pathogen detection with 30- compared to 20-milliliter blood culture draws. J Clin Microbiol 2011;49:4047-51.
115. Kaditis AG, O'Marcaigh AS, Rhodes KH, Weaver AL, Henry NK. Yield of positive blood cultures in pediatric oncology patients by a new method of blood culture collection. Pediatr Infect Dis J 1996;15:615-20.
116. Gonsalves WI, Cornish N, Moore M, Chen A, Varman M. Effects of volume and site of blood draw on blood culture results. J Clin Microbiol 2009;47:3482-5.
117. Kellogg JA, Manzella JP, Bankert DA. Frequency of low-level bacteremia in children from birth to fifteen years of age. J Clin Microbiol 2000;38:2181-5.
118. Mirrett S, Reller LB, Petti CA, Woods CW, Vazirani B, Sivadas R, et al. Controlled clinical comparison of BacT/ALERT standard aerobic medium with BACTEC standard aerobic medium for culturing blood. J Clin Microbiol 2003;41:2391-4.
119. Roh KH, Kim JY, Kim HN, Lee HJ, Sohn JW, Kim MJ, et al. Evaluation of BACTEC Plus aerobic and anaerobic blood culture bottles and BacT/Alert FAN aerobic and anaerobic blood culture bottles for the detection of bacteremia in ICU patients. Diagn Microbiol Infect Dis 2012;73:239-42.
120. Flayhart D, Borek AP, Wakefield T, Dick J, Carroll KC. Comparison of BACTEC PLUS blood culture media to BacT/Alert FA blood culture media for detection of bacterial pathogens in samples containing therapeutic levels of antibiotics. J Clin Microbiol 2007;45:816-21.
121. Miller NS, Rogan D, Orr BL, Whitney D. Comparison of BD Bactec Plus blood culture media to VersaTREK Redox blood culture media for detection of bacterial pathogens in simulated adult blood cultures containing therapeutic concentrations of antibiotics. J Clin Microbiol 2011;49:1624-7.
122. Ziegler R, Johnscher I, Martus P, Lenhardt D, Just HM. Controlled clinical laboratory comparison of two supplemented aerobic and anaerobic media used in automated blood culture systems to detect bloodstream infections. J Clin Microbiol 1998;36:657-61.
123. Wilson ML, Mirrett S, Meredith FT, Weinstein MP, Scotto V, Reller LB. Controlled clinical comparison of BACTEC plus anaerobic/F to standard anaerobic/F as the anaerobic companion bottle to plus aerobic/F medium for culturing blood from adults. J Clin Microbiol 2001;39:983-9.
124. Knepper JG, Anthony BF. Diminished growth of *Pseudomonas aeruginosa* in unvented blood-culture bottles. Lancet 1973;2:285-7.
125. Gantz NM, Medeiros AA, Swain JL, O'Brien TF. Vacuum blood-culture bottles inhibiting growth of *Candida* and fostering growth of *Bacteroides*. Lancet 1974;2:1174-6.

126. Mirrett S, Everts RJ, Reller LB. Controlled comparison of original vented aerobic fan medium with new nonvented BacT/ALERT FA medium for culturing blood. J Clin Microbiol 2001;39:2098-101.
127. Waites KB, Canupp KC. Evaluation of BacT/ALERT system for detection of *Mycoplasma hominis* in simulated blood cultures. J Clin Microbiol 2001;39:4328-31.
128. Eng J, Holten E. Gelatin neutralization of the inhibitory effect of sodium polyanethol sulfonate on *Neisseria meningitidis* in blood culture media. J Clin Microbiol 1977;6:1-3.
129. Davies S, Eggington R. Recovery of *Mycoplasma hominis* from blood culture media. Med Lab Sci 1991;48:110-3.
130. Courcol RJ, Durocher AV, Roussel-Delvallez M, Fruchart A, Martin GR. Routine evaluation of BACTEC NR-16A and NR-17A media. J Clin Microbiol 1988;26:1619-22.
131. Smith SM, Eng RH. Effectiveness of antibiotic removal by the antibiotic-binding blood culture systems. Diagn Microbiol Infect Dis 1985;3:201-12.
132. Wright AJ, Thompson RL, McLimans CA, Wilson WR, Washington JA. The antimicrobial removal device. A microbiological and clinical evaluation. Am J Clin Pathol 1982;78:173-7.
133. Jessamine PG, Hoban DJ, Forward KR. Positive Bactec resin cultures do not influence antimicrobial selection. Diagn Microbiol Infect Dis 1990;13:281-4.
134. Kim TJ, Weinstein MP. Update on blood cultures: how to obtain, process, report, and interpret. Clin Microbiol Infect 2013;19:513-20.
135. Washington JA. Collection, transport, and processing of blood cultures. Clin Lab Med 1994;14:59-68.
136. Hawkins BL, Peterson EM, de la Maza LM. Improvement of positive blood culture detection by agitation. Diagn Microbiol Infect Dis 1986;5:207-13.
137. Prag J, Nir M, Jensen J, Arpi M. Should aerobic blood cultures be shaken intermittently or continuously? APMIS 1991;99:1078-82.
138. Yagupsky P, Peled N, Press J, Abramson O, Abu-Rashid M. Comparison of BACTEC 9240 Peds Plus medium and isolator 1.5 microbial tube for detection of *Brucella melitensis* from blood cultures. J Clin Microbiol 1997;35:1382-4.
139. Herne V, Nelovkov A, Kutt M, Ivanova M. Diagnostic performance and therapeutic impact of LightCycler SeptiFast assay in patients with suspected sepsis. Eur J Microbiol Immunol (Bp) 2013;3:68-76.
140. Sullivan KV, Turner NN, Roundtree SS, McGowan KL. Rapid detection of methicillin-resistant *Staphylococcus aureus* (MRSA) and methicillin-susceptible *Staphylococcus aureus* (MSSA) using the KeyPath MRSA/MSSA blood culture test and the BacT/ALERT system in a pediatric population. Arch Pathol Lab Med 2013;137:1103-5.
141. Buchan BW, Ginocchio CC, Manii R, Cavagnolo R, Pancholi P, Swyers L, et al. Multiplex identification of gram-positive bacteria and resistance determinants directly from positive blood culture broths: evaluation of an automated microarray-based nucleic acid test. PLoS Med 2013;10:e1001478.
142. Martiny D, Debaugnies F, Gateff D, Gerard M, Aoun M, Martin C, et al. Impact of rapid microbial identification directly from positive blood cultures using matrix-assisted laser desorption/ionization time-of-flight mass spectrometry on patient management. Clin Microbiol Infect 2013;19:E568-E581.
143. Yossepowitch O, Dan M, Kutchinsky A, Gottesman T, Schwartz-Harari O. A cost-saving algorithm for rapid diagnosis of *Staphylococcus aureus* and susceptibility to oxacillin directly from positive blood

culture bottles by combined testing with BinaxNOW(R) *S. aureus* and Xpert MRSA/SA assay. *Diagn Microbiol Infect Dis* 2014.

144. Lagace-Wiens PR, Alfa MJ, Manickam K, Karlowsky JA. Thermostable DNase is superior to tube coagulase for direct detection of *Staphylococcus aureus* in positive blood cultures. *J Clin Microbiol* 2007;45:3478-9.
145. Qian Q, Eichelberger K, Kirby JE. Rapid identification of *Staphylococcus aureus* in blood cultures by use of the direct tube coagulase test. *J Clin Microbiol* 2007;45:2267-9.
146. Myrick BA, Ellner PD. Evaluation of the latex slide agglutination test for identification of *Staphylococcus aureus*. *J Clin Microbiol* 1982;15:275-7.
147. Davies S, Gear JE, Mason CM, McIntyre SM, Hall L. Streptococcus grouping latex kits: evaluation of five commercially available examples. *Br J Biomed Sci* 2003;60:136-40.
148. McDonald CL, Chapin K. Rapid identification of *Staphylococcus aureus* from blood culture bottles by a classic 2-hour tube coagulase test. *J Clin Microbiol* 1995;33:50-2.
149. Enoch DA, Cooke FJ, Guha S, Brown NM. Thermostable nuclease: a study of clinical impact. *J Antimicrob Chemother* 2008;61:754-5.
150. Speers DJ, Olma TR, Gilbert GL. Evaluation of four methods for rapid identification of *Staphylococcus aureus* from blood cultures. *J Clin Microbiol* 1998;36:1032-4.
151. Jin WY, Jang SJ, Lee MJ, Park G, Kim MJ, Kook JK, et al. Evaluation of VITEK 2, MicroScan, and Phoenix for identification of clinical isolates and reference strains. *Diagn Microbiol Infect Dis* 2011;70:442-7.
152. Mittman SA, Huard RC, Della-Latta P, Whittier S. Comparison of BD Phoenix to Vitek 2, Microscan MICroSTREP, and Etest for antimicrobial susceptibility testing of *Streptococcus pneumoniae*. *J Clin Microbiol* 2009;47:3557-61.
153. Carrara L, Navarro F, Turbau M, Seres M, Moran I, Quintana I, et al. Molecular diagnosis of bloodstream infections with a new dual priming oligonucleotide-based multiplex PCR assay. *J Med Microbiol* 2013.
154. Lyle N, Boyd J. The potential for PCR based testing to improve diagnosis and treatment of sepsis. *Curr Infect Dis Rep* 2013;15:372-9.
155. Carbonnelle E, Mesquita C, Bille E, Day N, Dauphin B, Beretti JL, et al. MALDI-TOF mass spectrometry tools for bacterial identification in clinical microbiology laboratory. *Clin Biochem* 2011;44:104-9.
156. Klein S, Zimmermann S, Kohler C, Mischnik A, Alle W, Bode KA. Integration of matrix-assisted laser desorption/ionization time-of-flight mass spectrometry in blood culture diagnostics: a fast and effective approach. *J Med Microbiol* 2012;61:323-31.
157. Wuppenhorst N, Consoir C, Lorch D, Schneider C. Direct identification of bacteria from charcoal-containing blood culture bottles using matrix-assisted laser desorption/ionisation time-of-flight mass spectrometry. *Eur J Clin Microbiol Infect Dis* 2012.
158. Dekker JP, Branda JA. MALDI-TOF mass spectrometry in the clinical microbiology laboratory. *Clinical Microbiology Newsletter* 2011;33:87-93.
159. Schmidt V, Jarosch A, Marz P, Sander C, Vacata V, Kalka-Moll W. Rapid identification of bacteria in positive blood culture by matrix-assisted laser desorption ionization time-of-flight mass spectrometry. *Eur J Clin Microbiol Infect Dis* 2012;31:311-7.

160. British Society for Antimicrobial Chemotherapy. BSAC Methods for Antimicrobial Susceptibility Testing. 2013.
161. Coyle MB, McGonagle LA, Plorde JJ, Clausen CR, Schoenknecht FD. Rapid antimicrobial susceptibility testing of isolates from blood cultures by direct inoculation and early reading of disk diffusion tests. J Clin Microbiol 1984;20:473-7.
162. Mushtaq S, Warner M, Cloke J, Afzal-Shah M, Livermore DM. Performance of the Oxoid M.I.C.Evaluator Strips compared with the Etest assay and BSAC agar dilution. J Antimicrob Chemother 2010;65:1702-11.
163. NHS Quality Improvement Scotland. Blood Culture: Driver Diagram, Implementation Framework, Priority Elements Checklist, Data Measurement Tool and Measurement Plan Working Draft. 2011.
164. Harvey DJ, Albert S. Standardized definition of contamination and evidence-based target necessary for high-quality blood culture contamination rate audit. J Hosp Infect 2013;83:265-6.
165. Freeman JT, Chen LF, Sexton DJ, Anderson DJ. Blood culture contamination with Enterococci and skin organisms: implications for surveillance definitions of primary bloodstream infections. Am J Infect Control 2011;39:436-8.
166. Roth A, Wiklund AE, Palsson AS, Melander EZ, Wullt M, Cronqvist J, et al. Reducing blood culture contamination by a simple informational intervention. J Clin Microbiol 2010;48:4552-8.
167. European Parliament. UK Standards for Microbiology Investigations (SMIs) use the term "CE marked leak proof container" to describe containers bearing the CE marking used for the collection and transport of clinical specimens. The requirements for specimen containers are given in the EU *in vitro* Diagnostic Medical Devices Directive (98/79/EC Annex 1 B 2.1) which states: "The design must allow easy handling and, where necessary, reduce as far as possible contamination of, and leakage from, the device during use and, in the case of specimen receptacles, the risk of contamination of the specimen. The manufacturing processes must be appropriate for these purposes".
168. Official Journal of the European Communities. Directive 98/79/EC of the European Parliament and of the Council of 27 October 1998 on *in vitro* diagnostic medical devices. 7-12-1998. p. 1-37.
169. Klaerner HG, Eschenbach U, Kamereck K, Lehn N, Wagner H, Miethke T. Failure of an automated blood culture system to detect nonfermentative gram-negative bacteria. J Clin Microbiol 2000;38:1036-41.
170. Seegmuller I, Eschenbach U, Kamereck K, Miethke T. Sensitivity of the BacT/ALERT FA-medium for detection of *Pseudomonas aeruginosa* in pre-incubated blood cultures and its temperature-dependence. J Med Microbiol 2004;53:869-74.
171. Narasimhan SL, Weinstein AJ. Infective endocarditis due to a nutritionally deficient streptococcus. J Pediatr 1980;96:61-2.
172. McCarthy LR, Bottone EJ. Bacteremia and endocarditis caused by satelliting streptococci. Am J Clin Pathol 1974;61:585-91.
173. Casetta A, Derouin V, Boussougant Y. Absence of spontaneous autolysis of *Streptococcus pneumoniae* in aerobic fan culture bottles in a commercial blood culture system. Eur J Clin Microbiol Infect Dis 1996;15:616-7.
174. Davis JC. A gram stain for smears of blood cultures, body fluids and tissues. Am J Med Technol 1976;42:417-23.

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175. Health and Safety Executive. Safe use of pneumatic air tube transport systems for pathology specimens. 9/99.
176. Department for transport. Transport of Infectious Substances, 2011 Revision 5. 2011.
177. World Health Organization. Guidance on regulations for the Transport of Infectious Substances 2013-2014. 2012.
178. Home Office. Anti-terrorism, Crime and Security Act. 2001 (as amended).
179. Advisory Committee on Dangerous Pathogens. The Approved List of Biological Agents. Health and Safety Executive. 2013. p. 1-32
180. Advisory Committee on Dangerous Pathogens. Infections at work: Controlling the risks. Her Majesty's Stationery Office. 2003.
181. Advisory Committee on Dangerous Pathogens. Biological agents: Managing the risks in laboratories and healthcare premises. Health and Safety Executive. 2005.
182. Advisory Committee on Dangerous Pathogens. Biological Agents: Managing the Risks in Laboratories and Healthcare Premises. Appendix 1.2 Transport of Infectious Substances - Revision. Health and Safety Executive. 2008.
183. Centers for Disease Control and Prevention. Guidelines for Safe Work Practices in Human and Animal Medical Diagnostic Laboratories. MMWR Surveill Summ 2012;61:1-102.
184. Health and Safety Executive. Control of Substances Hazardous to Health Regulations. The Control of Substances Hazardous to Health Regulations 2002. 5th ed. HSE Books; 2002.
185. Health and Safety Executive. Five Steps to Risk Assessment: A Step by Step Guide to a Safer and Healthier Workplace. HSE Books. 2002.
186. Health and Safety Executive. A Guide to Risk Assessment Requirements: Common Provisions in Health and Safety Law. HSE Books. 2002.
187. Health Services Advisory Committee. Safe Working and the Prevention of Infection in Clinical Laboratories and Similar Facilities. HSE Books. 2003.
188. British Standards Institution (BSI). BS EN12469 - Biotechnology - performance criteria for microbiological safety cabinets. 2000.
189. British Standards Institution (BSI). BS 5726:2005 - Microbiological safety cabinets. Information to be supplied by the purchaser and to the vendor and to the installer, and siting and use of cabinets. Recommendations and guidance. 24-3-2005. p. 1-14
190. Baron EJ, Miller JM, Weinstein MP, Richter SS, Gilligan PH, Thomson RB, Jr., et al. A Guide to Utilization of the Microbiology Laboratory for Diagnosis of Infectious Diseases: 2013 Recommendations by the Infectious Diseases Society of America (IDSA) and the American Society for Microbiology (ASM). Clin Infect Dis 2013;57:e22-e121.
191. Kelly MT, Roberts FJ, Henry D, Geere I, Smith JA. Clinical comparison of isolator and BACTEC 660 resin media for blood culture. J Clin Microbiol 1990;28:1925-7.
192. Bosshard PP. Incubation of fungal cultures: how long is long enough? Mycoses 2011;54:e539-e545.

193. Morris AJ, Byrne TC, Madden JF, Reller LB. Duration of incubation of fungal cultures. J Clin Microbiol 1996;34:1583-5.
194. Bhatti S, Vilenski L, Tight R, Smego RA, Jr. Histoplasma endocarditis: clinical and mycologic features and outcomes. J Infect 2005;51:2-9.
195. Public Health England. Laboratory Reporting to Public Health England: A Guide for Diagnostic Laboratories. 2013. p. 1-37.
196. Department of Health. Health Protection Legislation (England) Guidance. 2010. p. 1-112.
197. Lai CC, Wang CY, Liu WL, Huang YT, Liao CH, Hsueh PR. Time to positivity in blood cultures of staphylococci: clinical significance in bacteremia. J Infect 2011;62:249-51.